

The essence of nature in expert hands.

BEGO IMPLANT SYSTEMS

naturesQue ColTect P

Assuming that the clinical requirements for a barrier membrane are satisfied, individual preferences for the handling properties play a major role in the choice of collagen membrane. The easy handling of naturesQue *ColTect P* enables diverse applications in therapeutic concepts that include the use of regenerative measures.

The properties of naturesQue ColTect P

Structure, design, and stability

- Produced from porcine peritoneum
- Stabilization of collagen fibers using minimal crosslinking prevents premature enzymatic cleavage of the collagen fibers
- Reproducible barrier function (12–16 weeks)

Handling and feel

- Membrane can be applied wet or dry
- Feel when dry: Stable and rigid, sturdy
- Feel and properties when wetted: Can be draped, sturdy, adhesive, not sticky, can be repositioned
- Not side specific
- Can be fixed with pins

naturesQue ColTect P	
Item no.	Size
58496	15×20 mm
58497	20×30 mm
58498	30×40 mm
Manufacturer: Collagen Matrix Inc.	

